ADAPTARAP™

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ADAPTARAP™

Instructions for Use

Once it is determined that a wrap should be applied, the following steps will help achieve the best performance with ADAPTARAP.

If a wound treatment is to be applied, it is a best practice to use gauze to hold the treatment in place. This is critical if the treatment is a spray, ointment or other liquid which could prematurely dissolve ADAPTARAP.

Begin by applying ADAPTARAP in a similar manner as traditional footwraps. Apply with constant pressure while wrapping the foot, reducing gaps between layers to promote better adherence. The non-elastic nature ensures a close fit without increasing the risk of pressure wounds.

Avoid completing multiple 360 degree turns around the hoof without passing between the toes each time. This will accelerate the process of self-removal. ADAPTARAP should be applied over the heels and up the back of the leg to just below the dewclaws to prevent the wrap from sliding off during the treatment period.

ADAPTARAP was designed to remain on the hoof between two hours and five days of application, when properly applied. A minimum of two turns covering each of the heels would be needed to achieve a minimum persistence. Persistency can be extended by increasing the number of turns and length of ADAPTARAP that is applied.

After applying the desired amount of ADAPTARAP, wet the tail end of the wrap with a small amount of water and firmly press it against the rest of the wrap, creating a seal. A spray bottle works well to lightly wet the wrap. In times of extreme cold temperatures, use warm water for best results.